




**NEW
DESIGN**

DISPOSABLE URETERAL GUIDE SHEATH

Disposable Ureteral Guide Sheath



SCOPE OF APPLICATION

This product is used in endoscopic urological surgery to help dilate the ureter to build access for insertion and removal of the endoscope and associated surgical instruments as well as fluid injection.



PRODUCT FEATURES

The world's first cross-shaped negative pressure suction structure design

- This pioneering design is more in line with ergonomics, making the operation of negative pressure suction stone more convenient.

Detachable Design


- The sheath tube and the negative pressure connector are designed to be detachable. This pioneering design can adjust the angle of the negative pressure joint at will according to the human body position or operation habits, which brings higher convenience to the operation.

Hydrophilic coating

- Hydrophilic coating at both inner and outer tube to avoid damage to the urinary tract and facilitate the discharge of calculus fragments.

Passive bending

- Flexible-Y™ front end facilitates to passively bend via endoscope to observe the stone in narrow renal calyx and improve visual field.



Soft and smooth design

- Flexible tip and smooth transition of connection port to protect ureter and device from damage during access.

Model differentiation	Model specification				
Cross-shaped	ZSR-UAC-1035A	ZSR-UAC-1040A	ZSR-UAC-1045A	ZSR-UAC-1050A	ZSR-UAC-1055A
	ZSR-UAC-1135A	ZSR-UAC-1140A	ZSR-UAC-1145A	ZSR-UAC-1150A	ZSR-UAC-1155A
	ZSR-UAC-1235A	ZSR-UAC-1240A	ZSR-UAC-1245A	ZSR-UAC-1250A	ZSR-UAC-1255A
	ZSR-UAC-1335A	ZSR-UAC-1340A	ZSR-UAC-1345A	ZSR-UAC-1350A	ZSR-UAC-1355A
	ZSR-UAC-1435A	ZSR-UAC-1440A	ZSR-UAC-1445A	ZSR-UAC-1450A	ZSR-UAC-1455A



Dongguan ZSR Biomedical Technology Company Limited

Add : Unit 448, Qingfeng Road, Taihu Village, Sanzhong Village Committee, Qingxi Town , 523651 DONGGUAN City Guangdong Province P.R. China
 Fax : +86 769-87328976 Postcode : 523651
 Email:sales@zsrbiomedical.com

www.zsrbiomedical.com

