


DISPOSABLE  
URETERAL GUIDE SHEATH

# Disposable Ureteral Guide Sheath




## SCOPE OF APPLICATION

This product is used in endoscopic urological surgery to help dilate the ureter to build access for insertion and removal of the endoscope and associated surgical instruments as well as fluid injection.

## PRODUCT FEATURES

### Soft and smooth design

- Flexible tip and smooth transition of connection port to protect ureter and device from damage during access.



smooth transition of connection port

### Hydrophilic coating

- Hydrophilic coating at both inner and outer tube to avoid damage to the urinary tract and facilitate the discharge of calculus fragments.

### Flexible-Y™ type negative suction sheath

- Flexible-Y™ front end facilitates to passively bend via endoscope to observe the stone in narrow renal calyx and improve visual field.

Clear the stone during smashing to save the surgery time, meanwhile, improve the stone clearance rate.




Flexible Sheath

### Security

- Rigid sheath, effectively reduce intra-operative flushing pressure and improve safety.

### Multi-specification available

- Meet the different needs of clinical practice.



### Accuracy

- Radiography tip, clear and accurate positioning.



Radiography tip

Model differentiation	Item	ID	OD	Length	Model differentiation	Item	ID	OD	Length
Straight type	ZSR-UA-1035	10Fr	12Fr	35cm	Y type	ZSR-UAY-0935	9Fr	11Fr	35cm
	ZSR-UA-1045	10Fr	12Fr	45cm		ZSR-UAY-0945	9Fr	11Fr	45cm
	ZSR-UA-1135	11Fr	13Fr	35cm		ZSR-UAY-1035	10Fr	12Fr	35cm
	ZSR-UA-1145	11Fr	13Fr	45cm		ZSR-UAY-1045	10Fr	12Fr	45cm
	ZSR-UA-1235	12Fr	14Fr	35cm		ZSR-UAY-1135	11Fr	13Fr	35cm
	ZSR-UA-1245	12Fr	14Fr	45cm		ZSR-UAY-1145	11Fr	13Fr	45cm
	ZSR-UA-1335	13Fr	15Fr	35cm		ZSR-UAY-1235	12Fr	14Fr	35cm
	ZSR-UA-1345	13Fr	15Fr	45cm		ZSR-UAY-1245	12Fr	14Fr	45cm
	ZSR-UA-1435	14Fr	16Fr	35cm		ZSR-UAY-1335	13Fr	15Fr	35cm
	ZSR-UA-1445	14Fr	16Fr	45cm		ZSR-UAY-1345	13Fr	15Fr	45cm
W type	ZSR-UAW-1035	10Fr	12Fr	35cm	ZSR-UAY-1435	14Fr	16Fr	35cm	
	ZSR-UAW-1045	10Fr	12Fr	45cm	ZSR-UAY-1445	14Fr	16Fr	45cm	
	ZSR-UAW-1135	11Fr	13Fr	35cm					
	ZSR-UAW-1145	11Fr	13Fr	45cm					
	ZSR-UAW-1235	12Fr	14Fr	35cm					
	ZSR-UAW-1245	12Fr	14Fr	45cm					
	ZSR-UAW-1335	13Fr	15Fr	35cm					
	ZSR-UAW-1345	13Fr	15Fr	45cm					
ZSR-UAW-1435	14Fr	16Fr	35cm						
ZSR-UAW-1445	14Fr	16Fr	45cm						



### Dongguan ZSR Biomedical Technology Company Limited

Add : Unit 448, Qingfeng Road, Taihu Village, Sanzhong Village Committee, Qingxi Town , 523651 DONGGUAN City Guangdong Province P.R. China  
 Fax : +86 769-87328976 Postcode : 523651  
 Email:sales@zsrbiomedical.com

www.zsrbiomedical.com

